

# ***Life Diagnostics, Inc.***

610-431-7707 • Fax: 610-431-7818

Info@lifediagnostics.com • www.lifediagnostics.com

## ***Anti-Human Myoglobin Monoclonal***

Catalog No.	10080
Subtype	IgG1
Purity	Protein G purified
Protein	1-10 mg/ml
Buffer	10 mM Sodium Phosphate, 150 mM NaCl, pH 7.3
Storage	2-8°C
Immunogen	Highly purified human cardiac myoglobin
Pairing Information	This monoclonal can be used in an immunoassay for specific determinations of human myoglobin. When used as a capture monoclonal it pairs well with Affinity Purified Goat Anti-Myoglobin 10060 used as a conjugate antibody

P.O. Box 5205 • West Chester, PA 19380 • USA