

# ***Life Diagnostics, Inc.***

610-431-7707 • Fax: 610-431-7818 • [Info@lifediagnostics.com](mailto:Info@lifediagnostics.com)  
[www.lifediagnostics.com](http://www.lifediagnostics.com)

## PRODUCT DATA SHEET

**PRODUCT:** Rat Serum Control

**CATALOG #:** RSCON01

**LOT #:** N-A1309B

**FORMAT:** Lyophilized Wistar rat serum.  
Reconstitute by adding 200 µl of distilled or deionized water.

**DESCRIPTION:** Contains the indicated levels of serum proteins (see table below)\*

**STORAGE:** Store lyophilized control at -20°C or colder

Analyte	Concentration (mg/ml)
Alpha-1-acid Glycoprotein	0.242 +/- 0.036
Alpha-2-Macroglobulin	0.98 +/- 0.183
C-Reactive Protein	0.436 +/- 0.173
Ceruloplasmin	1.245 +/- 0.231
Clusterin	0.026 +/- 0.004
Haptoglobin	1.05 +/- 0.325
Hemopexin	1.38 +/- 0.35
Murinoglobulin	5.44 +/- 1.06
Thiostatin	0.966 +/- 0.098
IgG	6.47 +/- 1.12
IgM	5.20 +/- 0.47

P.O. Box 5205 • West Chester PA 19380 • USA